

# Dynamic Light Scattering


Shambo Sinha – 20CH10060

Utkarsh Baghla – 20CH10090


Sneha Majumder – 20CH10064

Priya Nandi – 20CH10046

Sambhuti – 20CH10054




# Introduction – Dynamic Light Scattering



# Introduction – Dynamic Light Scattering

.....Contd.


- ❖ The suspended particles of the colloidal dispersion under investigation undergo Brownian motion.
- ❖ This motion results in fluctuations in the distances between the particles and hence also in fluctuations of the phase relations of the scattered light.
- ❖ Due to the random motion of the suspended particles within the sample, the interference can be either constructive or destructive.
- ❖ The number of particles within the scattering volume may vary in time
- ❖ The net result is fluctuating intensity
- ❖ Scattered waves interference generates a net scattered light intensity  $I_s(t)$ .



# Dynamic Light Scattering – Time Correlation

---

- ❖ Spectral characteristics of the scattered lights depend on motion of scatters in time scale
- ❖ The objective of scattered light experiment is to measure time correlation functions of scattered field and/or scattered intensity
- ❖ Time correlation functions and their spectral densities are central to light scattering experiments
- ❖ Due to erratic thermal molecular motion, scattered field varies randomly at the detector.
- ❖ Time dependent correlation functions related to noise and fluctuation are relevant for the study of DLS





In general, autocorrelation function decays like exponential function

Where  $\tau_r$  is the relaxation time or correlation time

The autocorrelation function is the measure of similarities two signals, namely  $A(t)$  and  $A(t + \tau)$ .

# DLS – Particle Sizing Governing Equations

---

The corresponding measured normalized intensity correlation function is written as:

$$g_2(q, \tau) = \frac{\langle I_s(q, t) I_s(q, t + \tau) \rangle}{\langle |I_s(q, t)|^2 \rangle},$$

By means of the Siegert relation can be related to the electric field correlation function  $g_1(q, \tau)$ :

$$g_2(q, \tau) = 1 + \beta |g_1(q, \tau)|^2$$

$\beta$  being the so-called intercept and where the field correlation function  $g_1(q, \tau)$  is defined as

$$g_1(q, \tau) = \frac{\langle E_s(q, t) E_s^*(q, t + \tau) \rangle}{\langle |E_s(q, t)|^2 \rangle}$$

The field correlation function may be used to determine the diffusion coefficient  $D$  of the scatterers. For a monodisperse sample, is fitted to an exponential function

$$g_1(q, \tau) = \exp(-\Gamma \tau)$$

yielding the decay rate  $\Gamma$ . From its definition:


$$\Gamma = q^2 D$$

one obtains the diffusion coefficient  $D$  and by using the Stokes-Einstein equation one can calculate the hydrodynamic radius

$$R = \frac{kT}{6\pi\eta D}$$

# DLS – Particle Sizing Example

---



Source: Jia, Z.; Li, J.; Gao, L.; Yang, D.; Kanaev, A. Dynamic Light Scattering: A Powerful Tool for In Situ Nanoparticle Sizing. *Colloids Interfaces* **2023**, 7, 15.  
<https://doi.org/10.3390/colloids7010015>

# Particle Size Characterization in Colloids

---

**Colloids:** Colloids are mixtures in which one substance is dispersed as relatively large solid particles or liquid droplets throughout a solid, liquid, or gaseous medium.

Colloids are crucial in paints, pigments, pharmaceuticals, food, cosmetics, ceramics, and personal care. Particle size, polydispersity, shape, and surface characteristics significantly influence formulation efficacy.

## Methods for measuring particle size:

### Microscopic Imaging

High resolution but requires counting numerous particles

### Electrical sensing (Coulter) counters

Quick, direct; challenges with non-spherical particles

### Hydrodynamic or field flow fractionation

Efficient, struggles with polydisperse samples

### Disc centrifuge particle sizing

High res, time-consuming and complex process

### Size exclusion chromatography

Suitable for macromolecules

### Scattering techniques

Valuable info, may require complex data analysis



# Scattering Techniques

## Static Light Scattering (SLS)

- Measures scattering intensity at various angles.
- Limited to particles similar or larger than the wavelength of light..

## Dynamic Light Scattering (DLS)


- Utilizes Brownian motion-induced intensity fluctuations.
- Measures dynamic size changes in particles undergoing random motion.
- Ensemble-averaged submicron size estimates.
- Fast data acquisition

## Analyzing Particle Size via Intensity Fluctuations

- Utilizes speckle patterns resulting from scattered light
- Measures temporal fluctuations in intensity to determine particle size
- Autocorrelation function helps determine particle diffusion coefficients
- Characterizes intensity fluctuations by comparing signal at different times

**Autocorrelation Function:** 
$$g^2(\tau) = \frac{\langle I(t)I(t + \tau) \rangle}{\langle I(t) \rangle^2}$$

When  $\tau$  equals zero,  $g^2(\tau)$  equals 1, indicating that the signal is perfectly correlated at this point, reflecting uniform particle movement




# ANALYSIS OF AUTOCORRELATION FUNCTION

## Objective

- Compare the particle sizes based on an analysis of the graph depicting the correlation as a function of delay time

## Procedure (Simplified)


- The data set for the large and small particles has been extracted from the given graph.
- The dataset was subjected to a log transformation and then plotted to calculate the slope of the graph.




# ANALYSIS OF AUTOCORRELATION FUNCTION

.....Contd.

$\ln(G')$  vs Time - smaller particles



$\ln(G')$  vs Time - larger particles



# ANALYSIS OF AUTOCORRELATION FUNCTION

.....Contd.

For a suspension of spherical particles undergoing Brownian diffusion, the autocorrelation function decays exponentially with the delay time  $\tau$  and is given as:

$$g^1(\tau) = A \cdot e^{-Dq^2\tau} + B$$

where A: Amplitude of the correlation function

B: Baseline, accounts for any background or noise in the measurements

D: Translational diffusion coefficient of the particles, quantifies the ability of the particles to diffuse or spread out in the suspension. A higher D value implies faster diffusion of the particles.

q is the magnitude of the scattering vector given as  $(4\pi n/\lambda) \sin(\theta/2)$ . Here, n is the refractive index of the medium.

**Slope for smaller particles: -0.2990278052**

**Slope for larger particles: -0.09175909456**

For spherical particles, the hydrodynamic radius  $R_h$  can be obtained from the translational diffusion coefficient ( $D$ ) using the Stokes–Einstein relationship:

$$D = kT / (6\pi\eta R_h)$$

where  $k$  is the Boltzmann's constant,  $\eta$  is the solvent viscosity, and  $T$  is the absolute temperature

From the slopes, assuming the values of  $T$ , viscosity and  $k$ , we get the values of  $R$  as:

$R$ (smaller particles) : 3.344170618

$R$  (larger particles) : 10.89810231

*Thank You*

